Tutorial 12: Centrifuge

Gavin Hearne Mohamm<u>ad Saleh Refahi</u>



Background

What is Centrifuge?

- Centrifuge is a metagenomic abundance estimation tool developed by researchers at JHU that labels and quantifies metagenomic reads
- Most metagenomic classifiers either suffer from either speed (such as Naive Bayes Classifier, PhymmBL, and MegaBLAST), or index size (kraken) issues. This tool aims to solve both of these
- This is done through the use of a novel indexing scheme based on two data structures:
 - Burrows-Wheeler transform
 - o Ferragina-Manzini (FM) index


Database sequence compression

Multiple genomes of the same species are compressed by storing near-identical sequences only once.

- This can reduce total sequence lengths by up to 89%.

An FM-index is then created based on these precompressed sequences

- The FM-index is quite small when compared to k-mer indexing methods that store all k-mers (such as kraken)




abaaba


The FM-index

The FM index combines the BTW with a few auxiliary data structures - this consists of 4 components:

- First column (F) of the BW matrix
- Last Column (L) of the BW matrix
- Checkpoints
 - Allows the steps of FM querying to be performed in O(1) time
- Suffix Array
 - Quickly finds where all occurrences of a query appear in the initial sequence

This has the potential to significantly reduce the combined size of the stored genomes, while also allowing for faster queries of the index





Building Indexes

Prebuilt:

The researchers behind Centrifuge have put together three main indexes that are readily available:

- p+h+v: bacterial, human, and viral genomes [~12G]
- p_compressed: bacterial genomes
 compressed at the species level [~4.2G]
- p_compressed+h+v: combination of the two above [~8G]

Custom:

If these indexes do not suit your needs, it is possible to build custom indexes from arbitrary sequences.

For each sequence, centrifuge needs the following:

- nodes.dmp file from the NCBI taxa dump to build the taxa tree
- sequence ID to taxonomy ID map



The FM-index allows for Centrifuge to exploit the advantages of both large and small *k*-mers:


- As *k*-mers size increases, precision is increased but sensitivity is reduced


Centrifuge begins by finding 16-bp exact matches, then extending those matches as far as possible. This is performed both with the forward and reverse complement.

Classifications are then performed using only mappings with at least one 22-bp match.

Higher taxonomic level classifications are performed by traversing up the tree and counting number of matches per rank







Abundance analysis

Centrifuge also is capable of performing abundance analysis at any taxonomic rank

- Likelihood for a specific configuration of species abundance a given the read assignments

To find the abundance that maximizes the likelihood function, the following EM procedure is repeated until the difference between the previous estimate of abundances and the current estimate is less than 10–10

- Expectation (E-step): the estimated number of reads assigned to species j.
- Maximization (M-Step): the updated estimate of species j's abundance

$$L(\alpha|C) = \prod_{i=1}^{R} \sum_{j=1}^{S} \frac{\alpha_{j} l_{j}}{\sum_{k}^{S} \alpha_{k} l_{k}} C_{ij},$$

$$n_j = \sum_{i=1}^R \frac{\alpha_j C_{ij}}{\sum_{k=1}^S \alpha_k C_{ik}}$$

$$\alpha_j' = \frac{n_j/l_j}{\sum_{k=1}^S n_k/l_k},$$

Usage

Installing Requirements Packages


```
one-time set up of Bioconda with the following commands.
[ ]: !conda config --add channels defaults
     !conda config --add channels bioconda
     !conda config --add channels conda-forge
     !conda config --set channel_priority strict
     Install the following packages in conda environments
[ ]: # Install mamba in the base conda environment
     # Note: This step might be optional, as mamba is not always necessary
     !conda install -n base --override-channels -c conda-forge mamba 'python abi=*=*cp*'
[ ]: # Install the ete3 package using conda
     !conda install ete3 --quiet --yes
     !conda install -c conda-forge biopython
     # Install the centrifuge package from the bioconda channel using conda
     !conda install bioconda::centrifuge
     # Install r-remotes and bioconductor-rsamtools using mamba
     !mamba install --yes --quiet r-remotes bioconductor-rsamtools
     # Install ipywidgets using pip
     !pip install ipywidgets
```

!pip install epi2melabs



Downloading Centrifuge Index

h+p+v+c: human genome, prokaryotic genomes, and viral genomes including SARS-CoV-2 genomes.



Centrifuge Command Parameters

Index Selection

-x <cf-idx>: Index filename prefix (minus trailing .X.cf)

Input Options

-1 <m1> -2 <m2>: Paired-end input files

-U <r>: Unpaired reads

Output Configuration

-S <filename>: File for classification output (stdout by default)

--report-file <report>: Tabular report output file (centrifuge_report.tsv by default)

Classification Settings

--host-taxids <taxids>: Preferred taxonomic IDs in classification

--exclude-taxids <taxids>: Excluded taxonomic IDs in classification

Performance


-p/--threads <int>: Number of alignment threads (default: 1)

--mm: Use memory-mapped I/O for index; suitable for multiple instances

Pavian: Interactive Analysis of Metagenomics Data



Pavian is a interactive browser application for analyzing and visualization metagenomics classification results from classifiers such as Kraken, Kraken_Uniq, Kraken 2, Centrifuge and MetaPhIAn



Pavian: Interactive Analysis of Metagenomics Data

```
!centrifuge-kreport -x hpvc results.txt > read_classifications.tsv.kraken
                                                                                                                                                                                  □ ↑ ↓ 古 早 盲
Loading taxonomy ...
Loading names file ...
/home/sr3622/miniconda3/bin/centrifuge-inspect:24: DeprecationWarning: the imp module is deprecated in favour of importlib and slated for removal in Python 3.12; see the module's documentation for alternative
 import imp
/home/sr3622/miniconda3/bin/centrifuge-inspect:24: DeprecationWarning: the imp module is deprecated in favour of importlib and slated for removal in Python 3.12; see the module's documentation for alternative
 import imp
import ipywidgets as widgets
from epi2melabs.notebook import InputForm, InputSpec
pavian form = InputForm(
   InputSpec('port', 'Aux. EPI2ME Labs port', widgets.IntText(8889)))
pavian_form.display()
VBox(children=(HBox(children=(Label(value='Aux. EPI2ME Labs port', layout=Layout(width='150px')), interactive(...
# running Pavian web server
# !echo "Checkina pavian install..."
script = """
                                                                                                          Open Url: localhost:8889
ncpus=parallel::detectCores()
options(Ncpus=ncpus)
remotes::install_github("fbreitwieser/pavian", upgrade=T, quiet=T)"""
_script = os.path.expanduser("~/.pavian_install.R")
with open(_script, "w") as fh:
   fh.write(script)
!Rscript $_script
!echo "Done."
!echo "Running pavian..."
port = pavian form.port
!R -e "pavian::runApp(host='0.0.0.0', port="$port")"
Running pavian...
```

Centrifuge Output

Classification Output

```
readID seqID
                   taxID
                                  2ndBestScore
                                                  hitLength
                                                              queryLength numMatches
                           score
   NS500207:12:H04WYAFXX:3:21408:22104:9367
                                              NZ CP014768.1
                                                              1813821 81 81 24 27 2
                                              species 562 81
   NS500207:12:H04WYAFXX:3:21408:22104:9367
                                                             81
                                                                 24
                                                                     27
   NS500207:12:H04WYAFXX:4:21410:2455:20204
                                              unclassified
                                                                 0
                                                                     0
                                                                         0
                                                                             16 1
                                              unclassified
   NS500207:12:H04WYAFXX:2:21105:21680:16668
                                                                             21 1
                                              unclassified
                                                                             15 1
234 NS500207:12:H04WYAFXX:3:11609:20259:9826
```


Report output

```
taxID
                 taxRank genomeSize numReads
                                                 numUniqueReads
                                                                abundance
  name
 Cellvibrio 10
                         4576573 1
                                         0.0
                 genus
 Pseudomonas fluorescens 294 species 6604895 2
                                                     0.0
                         301 species 4686340 1
 Pseudomonas oleovorans
                                                     0.0
 Xanthomonas 338 genus
                                         0
                         23419814
                                             0.0
6 Neisseria
             482 genus
                         2223758 1
                                         0.0
```

Pavian: Sankey visualization

Let see the HTML output!

read_classifications.tsv.kraken



Resources

https://labs.epi2me.io/notebooks/Metagenomic_classification_tutorial.html?

Kim, D., Song, L., Breitwieser, F. P., & Salzberg, S. L. (2016). Centrifuge: rapid and sensitive classification of metagenomic sequences. *Genome research*, *26*(12), 1721-1729.

Breitwieser, F. P., & Salzberg, S. L. (2020). Pavian: interactive analysis of metagenomics data for microbiome studies and pathogen identification. Bioinformatics, 36(4), 1303-1304.